



Ambra Mobile App - iOS

1.



In the iOS App Store, search for “Ambra” and download the Ambra Health app.


2.



SET UP AMBRA HEALTH APP

1. Navigate to your **Settings**
2. Scroll down and click “Ambra”
3. Type in “pediatrix.ambrahealth.com” on the Host/Vanity field
4. Toggle on (enable) **Show Account Picker**
5. Select “**Highest Quality**” in the Image Quality viewer settings

3.




SIGNING IN

Navigate to and click the Ambra Health app

Scroll down to the “Patients” row and click “Continue”



You will be directed to select “PEDIATRIX” from a dropdown
Select **PEDIATRIX**




Enter your login credentials for accessing the PEDIATRIX Imaging Portal

Note: You will be prompted to enter the 2-factor authentication

Ambra Mobile App - iOS


4.



ACCESS IMAGING

Once you are logged in, you will be taken to the studies worklist

Select the **imaging** you would like to view



Once selected, your imaging will appear


IMAGE NAVIGATION

Go back to the worklist

Show or hide overlay

Adjust brightness of the image

Scroll between different images



App guidance

Reset to original image

Toggle cine mode of multi-frame series

Navigate to another series